

SmartJet™

Liquid and Spray Delivery Systems

SmartJet eliminates the guesswork of syringe mixing.



Predictability

- Precision controlled Ten-to-One fluid delivery
- Assures reproducible results

Control

- Sleek, pistol grip handle for comfort and control
- Single-handed operation for easy manipulation



SmartJet, Naturally

- Dual lumen design to prevent clogging at the delivery tip
- Packaged sterile and assembles in seconds
- More economical than conventional applicators



HARVEST®

Developing technologies for **accelerating healing, naturally™**



SmartJet™ Liquid and Spray Delivery Systems

- Fluid Delivered:**
- **Precisely and Reliably**
 - **Faster and Easier**
 - **Economically**

Applicator Components



Ordering Information

Product Code No.	Product Description
LK/2	 SmartJet Liquid Applicator Kit Sterile, single-use liquid applicator kit, includes: one (1) dual lumen 3.5inch (9 cm) tip, one (1) 1 cc syringe, one (1) 10 cc syringe, one (1) handle, and one (1) clip. 5/case.
SK/S	 SmartJet Spray Applicator Kit Sterile, single-use spray applicator kit, includes: one (1) spray tip, one (1) 1 cc syringe, one (1) 10 cc syringe, one (1) handle, and one (1) clip. 5/case.

**To arrange an evaluation or for more information,
call 877-8-HARVEST (toll free) or visit us at www.harvesttech.com**

HARVEST®

Developing technologies for accelerating healing, naturally™

Harvest Technologies Corp
40 Grissom Road, Plymouth, MA 02360
Phone (508)732-7500 Fax (508)732-0400

Harvest Technologies GmbH
Zehnfeldstr. 240a, D-81825 Munich, Germany
Phone +49(0)89 437778-0 Fax +49(0)89 437778-10

CE
0086
BSI Registered Firm
ISO 9001 Registered